

- Hospital grade power cords fit a variety of common infusion pumps
- Items are shipped from stock and include a 90-day warranty

Power Cords

Right Angle Cord



Left Angle Cord



Cord Length (ft)	Gauge	Plug Angle	AIV Part Number
8	14	Straight	PC1079401
8	14	Right	PC1079501
8	14	Left	PC1079601
18	14	Straight	PC1079402
25	14	Left	PC1079603
8	16	Straight	PCM8
8	16	Right	PC1079506
8	16	Left	PC1079606
18	16	Straight	PC1079405
18	16	Left	PC1079604
25	16	Straight	PC1079404
For Use with Baxter® Colleague® Pumps			
10	16	Left	PC10761
For Use with Alaris® Signature® Series and Medley® Pumps			
10	16	Straight	PC10775
For Use with GE® Dash Monitors			
6	16	Down	PC11724

Specialized Repair Services
Accessories
Parts

- Fetal Monitoring
- Pulse Oximetry
- Cardiac Monitoring
- Infusion Therapy

AIV, Inc.
7485 Shipley Avenue
Harmans, MD 21077 USA
p 410.787.1300
f 410.787.1337

**Call us today 800.990.2911 or
visit us at www.aiv-inc.com**

*The manufacturers listed are the holders of their respective names and/or trademarks, and are not to be taken as an endorsement or affiliation with AIV, Inc.